Articles

Cat. No

PL(Product Liability) Insurance

# DAIHAN ThermoStable<sup>™</sup> IS-10 Precise Shaking Incubators, Top Door-type, Orbital Motion with Universal Platform, Fuzzy Control, with or without illuminators & Recorder, up to 60°C, ±0.2°C, 30~250 rpm



Precise Shaking Incubator, Top Door Type "ThermoStable IS-10R"

- with Certificate & Traceability : Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Useful for Biological Cultures under various Temperature
- Patented Jog-Dial Control System
- CE Certified
- \* Precise Temperature Control by Digital Fuzzy Controller, up to 60°C
- \* Digital LCD with Back-Light Function
- \* Minimum Noise and No Vibration
- \* Advanced Shaking Mechanism Provides Quiet Shaking and Precise Speed Control
- \* Automatic Stop of Shaking, Heating and Refrigeration when a Door opens
- \* Universal Spring Rack for Various Shapes/Sizes of Flasks (Optional)
- \* RS232C Interface for Monitoring and Controlling with PC

#### Specification

Model		ThermoStable IS-10	ThermoStable IS-10R	ThermoStable IS-10RL
Illumination System		None		FPL 55W×4ea
Temp. Range		Ambient +5℃ ~ 60℃	10°C ~ 60°C	
Temp. Accuracy & Uniformity		Accuracy : $\pm 0.2^{\circ}$ at 37°C, Uniformity : $\pm 0.5^{\circ}$ at 37°C		
Temp. Sensor		PT100 Sensor		
Heater		800 W	1 kW	1.2 kW
Compressor		None	1/3 HP	
Compressor Control		None	Intelligent Control of Refrigerator for preventing Overload of Compressor	
Refrigerant		CFC-Free (R-404A) Refrigeration System		
Shaking Speed & Stroke		30 ~ 250 rpm, 25mm Orbital Motion		
Timer & Alarm		99hr 59min (with Continuous run.), Error status & Timer-end		
Display		Digital LCD with Back-Light Function		
Controller		Digital Fuzzy Controller with Jog-Dial Switch (Turn+Push)		
Room Lamp		Krypton Lamp, 60W		None
Platform		Included Universal Platform for Various Shapes & Sizes of Flasks		
Optional Flask Holder		Durable & Flexible Plastic Flask Holder (up to 1000ml)		
Available Max. for	50~100ml	Available 68 ea Flask Holder "FH100", for 50~100ml Flask		
	200~300ml	Available 36 ea Flask Holder "FH250", for 200~300ml Flask		
Flask Holder for	500~1,000ml	Available 21 ea Flask Holder "FH500", for 500~1,000ml Flask		
Safety Circuit		Over Temp./ Current Protector, Leakage Breaker, Sensor Error Detector		
Material		Internal : Stainless steel (#304), External : Powder Coated steel		
Dimension		Internal : $750 \times 720 \times h400$ mm, External : $1,065 \times 900 \times h980$ mm		
Packing Size (w×d×h) (mm) & Gross Weight		1,220×1,060×h1,200mm		
		231 kg	238 kg	238 kg
Power Consumption		800 W	1.8 kW	2.1 kW
Power Supply & Cord/Plug*		1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz		

\* Other Specifications are available upon Customer's Request.

#### Ordering Information

#### Ambient +5°C ~60°C, With Universal Platform

DH.WIS00210Shaking Incubator, Top Door, 230V, "ThermoStable IS-10"DH.WIS00110, 120V, "ThermoStable IS-10"

#### $10^{\circ}C \sim 60^{\circ}C$ , With Universal Platform

DH.WIS04210Shaking Incubator, Top Door, 230V, "ThermoStable IS-10R"DH.WIS04110, 120V, "ThermoStable IS-10R"

#### $10^{\circ}C \sim 60^{\circ}C$ , with Universal Platform & Illumination System

DH.WIS05210 Shaking Incubator, Top Door, 230V, "ThermoStable IS-10RL" DH.WIS05110 , 120V, "ThermoStable IS-10RL"

### <Continued on next page>


# <Continued from... DAIHAN ThermoStable<sup>™</sup> IS-10 Precise Shaking Incubators >

## Optional Platform/Spring Rack/Flask Holders

DH.WIS504510	(1) Universal Spring Rack, w693 × d653 × h135mm, "SR510", for "IS-10(R/RL)"
DH.WIS504610	(2) Universal Platform, w688×h648mm, "SP610", for "IS-10(R/RL)"
DH.WIS504030	(3) Flask Holder, for 50~100ml Flask, "FH100", for "SHO", "SHR" and "IS"
DH.WIS504035	(4) Flask Holder, for 200~300ml Flask, "FH250", for "SHO", "SHR" and "IS"
DH.WIS504040	(5) Flask Holder, for 500~1000ml Flask, "FH500", for "SHO", "SHR" and "IS"
DH.WTH05100	(6) Temperature Recorder, with 6ea Ribbon Cassette & Chart Paper LED, 6ch/10sec, 144 × 144 × h233mm, "WTH5105"
DH.WTH05108 DH.WTH05105	Ribbon Cassette, 6 Colors, " <b>WTH5108</b> " Chart Paper, w100mm×L10m, 33 Days, " <b>WTH5105</b> "
	Validation (IQ, OQ) Service

DH.WIS31010	Validation Service(IQ, OQ), "ISV1010", for "ThermoStable IS-10"
DH.WIS31011	, "ISV1011", for "ThermoStable IS-10R"
DH.WIS31012	, "ISV1012", for "ThermoStable IS-10RL"





(1) Universal Spring Rack "SR510" (Optional)

(2) Universal Platform "SP610" (Included)



Flask Holders (Optional)



(6) Temperature Recorder : with 6ea Ribbon Cassette and Chart Paper

# WiseRemote® Remote Control Software

Connection of Wisd.22 Temperature Devices with Interface RS232C to a PC Export Recorded Data (Excel Spreadsheet)



RS232C Cable & Software



Precise Shaking Incubator, Top Door Type "ThermoStable IS-10R"